

IDC™ Soft Embolization Coil

PRECISION IN PLACEMENT

Soft Coil Design for Smooth & Elegant Finishing:


Thinner platinum wire and softer primary coil wind designed for smooth, tight nesting of the finishing coil during vessel embolization

Precise Placement:

Detachable, soft coil technology with interlocking arms allows for precise placement of the coil tail and reduced risk of catheter kickback

Clinical Versatility:

Range of sizes (down to 2 mm x 2 cm) facilitates use in extra small peripheral vessels and for delicate finishing needs



IDC™ Soft Embolization Coil

UPN	Description	Diameter (mm)	Length (cm)	Shape
M0033612020	IDC-18 Soft	2	2	2D
M0013612040	IDC-18 Soft	2	4	2D
M0013613060	IDC-18 Soft	3	6	2D
M0013613100	IDC-18 Soft	3	10	2D
M0013614040	IDC-18 Soft	4	4	2D
M0013614080	IDC-18 Soft	4	8	2D
M0013614120	IDC-18 Soft	4	12	2D
M0013615080	IDC-18 Soft	5	8	2D
M0013615120	IDC-18 Soft	5	12	2D

IDC Soft Embolization Coil – Available Sizes

		Coil Diameter			
		2 mm	3 mm	4 mm	5 mm
Coil Length	2 cm				
	4 cm				
	6 cm				
	8 cm				
	10 cm				
	12 cm				



IDC INTERLOCKING DETACHABLE COIL

CAUTION: Federal law (USA) restricts this device to sale by or on the order of a physician. Rx only. Prior to use, please see the complete "Directions for Use" for more information on Indications, Contraindications, Warnings, Precautions, Adverse Events, and Operator's Instructions.

INTENDED USE/INDICATIONS FOR USE: The IDC Coil is a modified interlocking detachable coil indicated to obstruct or reduce rate of blood flow in the peripheral vasculature. This device is not intended for neurovascular use. **CONTRAINDICATIONS:** None known. **GENERAL PRECAUTIONS:** Do not advance the IDC Interlocking Detachable Coil if it becomes lodged within the microcatheter. Determine the cause of the resistance and replace the microcatheter and coil if necessary. **ADVERSE EVENTS:** The complications that may result from a peripheral embolization procedure include, but are not limited to: • Complications related to catheterization (e.g., hematoma at the site of entry, vessel injury, etc.) • Death • Emboli • Foreign body reactions necessitating medical intervention • Hemorrhage • Infection necessitating medical intervention • Ischemia • Pain • Recanalization • Temporary neurological deficit • Tissue necrosis • Undesirable clot formation of the vasculature • Vasospasm **REV AA**

IDC is a registered or unregistered trademark of Boston Scientific Corporation or its affiliates.

Peripheral Interventions

300 Boston Scientific Way
Marlborough, MA 01752-1234
www.bostonscientific.com

To order product or for more information contact customer service at 1.888.272.1001.

© 2015 Boston Scientific Corporation or its affiliates. All rights reserved.

PI-299412-AA FEB2015