

# AZUR HydroPack™

Peripheral Coil System

## PACK WITHOUT COMPROMISE

In a category of its own



- Hydrogel expansion creates a solid core for mechanical occlusion<sup>1</sup>
- AZUR HydroPack™ coil can be deployed through 2.4Fr or 2.8Fr Microcatheters<sup>2</sup>
- Available in long lengths, up to 60cm<sup>1</sup>
- Optimized for trackability and microcatheter stability during coil delivery<sup>1</sup>

**TERUMO**  
INTERVENTIONAL  
SYSTEMS

# PUSHING BOUNDARIES

Terumo Interventional Systems is **committed to your success** with innovative procedural solutions and ongoing support for your most challenging cases.

We are relentlessly seeking new ways to help you achieve **better outcomes for more patients**.



## ORDERING INFORMATION

AZUR HYDROPACK™ Peripheral Coil System						
PRODUCT CODE	DESCRIPTION	LENGTH (CM)	MICROCATHETER COMPATIBILITY*		SOFTNESS	HYDROGEL
			PROGREAT 2.4Fr	PROGREAT 2.8Fr		
45-880005	AZUR HydroPack 18 5 cm	5	✓	✓	✓	✓
45-880010	AZUR HydroPack 18 10 cm	10	✓	✓	✓	✓
45-880020	AZUR HydroPack 18 20 cm	20	✓	✓	✓	✓
45-880035	AZUR HydroPack 18 35 cm	35	✓	✓	✓	✓
45-880050	AZUR HydroPack 18 50 cm	50	✓	✓	✓	✓
45-880060	AZUR HydroPack 18 60 cm	60	✓	✓	✓	✓

\*AZUR HydroPack is compatible with microcatheters ranging from 0.021" - 0.027" I.D.





AZUR HydroPack is designed to be a soft coil that finds and fills empty space<sup>1</sup>



Hydrogel expansion creates a solid core for mechanical occlusion<sup>1</sup>



AZUR HydroPack has a wider range of microcatheter compatibility vs Pod Packing Standard Coil System<sup>1</sup>

**FIND OUT MORE**  Phone: 800.888.3786  [terumo.com](http://terumo.com)

**Indications for Use:** The AZUR system is intended to reduce or block the rate of blood flow in vessels of the peripheral vasculature. It is intended for use in the interventional radiologic management of arteriovenous malformations, arteriovenous fistulae, aneurysms, and other lesions of the peripheral vasculature.

Manufacturer: MicroVention, Inc., 35 Enterprise, Aliso Viejo, CA 92656 USA  
Distributed by Terumo Medical Corporation, 265 Davidson Avenue, Somerset, NJ 08873.

**RX ONLY. Refer to the product labels and package insert for complete warnings, precautions, potential complications, and instructions for use.**

**References:**

1. Data on file.
2. AZUR HydroPack IFU100206

©2023 Terumo Medical Corporation. All brand names are trademarks or registered trademarks of their respective owners. PM-07185

